



tenda®

Bleed-No-More Paste™

Pulmonary Support Paste

Directions: Administer 1-60 cc syringe 2 to 4 hours prior to event. May also administer 1-60 cc syringe 6 hours out and another 2 to 4 hours out. Minimize horse's water intake until race is complete.

Ingredients: Honey, Apple Cider Vinegar, Dextrose, Glycerine, Xanthan Gum, Pantothenic Acid, Yunnan Baiyao, Corn Oil.

Precautions: Keep out of reach of children. For equine use only. Do not use on horses intended for food consumption. Do not mix with feed. Do not freeze. If accidentally frozen allow to thaw at room temperature for 24 hours. Not for use on pregnant mares.

Net Content: 60 cc

Manufactured in
the U.S.A by:

Tenda Horse Products
Sandusky, Ohio 44870
www.tendahorse.com

